

BIOMÉRIEUX

ENDONEXT™

Control Standard Endotoxin.



Your Ally in Advancing Quality

PIONEERING DIAGNOSTICS



CONTROL STANDARD ENDOTOXIN

ENDONEXT™ Standard is intended for use as a Bacterial Endotoxin Test (BET) reference reagent containing a defined amount of endotoxin from *Escherichia coli* O113:H10. Thanks to the final concentration of 10,000 EU/ml, it can be used as a quantitative quality control for finished pharmaceutical products, in-process control, research samples, and medical device testing.

- Ideal for hard spiking and recovery studies
- Made from the same strain as reference standard endotoxin
- Highly characterized batch-to-batch
- Specifically designed for ENDONEXT™ rFC assay portfolio



ORDERING INFORMATION

Reference	Description	Packaging
423583	ENDONEXT™ STANDARD (10,000 EU/ML)	6 bottles

Experience the quality results you deserve with the ENDONEXT™ offer. Learn more at www.biomerieux.com